

Marble Launch Lab**Day 1: Find the initial velocity**

Purpose: to calculate the initial velocity of a marble launched from a CPO the marble launcher at two different settings.

Procedures:

- 1) Materials: eye protection, CPO marble launcher, marble, measuring tape.
- 2) WEAR EYE PROTECTION
- 3) Set up your marble launcher by attaching it to a stand and setting the angle to 0° .
- 4) Measure the height of the launcher from the ground to the center of the launching tube.
- 5) Choose which two settings you will be measuring. The setting closest to the nozzle is considered setting 1. The setting closest to the back is considered setting 5.
- 6) At each setting launch the marble at least 4 times and measure the horizontal displacement of each trial.

Data:

Height: _____

Setting 1: _____	displacement
Trial 1	
Trial 2	
Trial 3	
Trial 4	
Average displacement	

Setting 2: _____	displacement
Trial 1	
Trial 2	
Trial 3	
Trial 4	
Average displacement	

Calculations:

1) Prove that $v_0 = \sqrt{\frac{\Delta x}{\frac{2h}{4.9}}}$

2) Find the initial velocity of the marble at each setting.

Results:

The initial velocity at setting _____ is _____.

The initial velocity at setting _____ is _____.

Day 2: Make Predictions

Purpose: To use what we know about projectile motion to predict where the marble is going to land when it is launched at two different angles.

Procedures:

- 1) Decide what two angles you are going to launch your marble at.
- 2) Calculate the displacement of the marble at each of the 4 setting combinations. You will have 4 different calculations to make: setting 1-angle 1, setting 1-angle 2, setting 2-angle 1, setting 2-angle 2.
- 3) Split the initial velocity into its components using $v_x = v_i \cos \theta$ and $v_y = v_i \sin \theta$
- 4) Calculate the time it will take to fall in y using $\Delta y = v_y t + \frac{1}{2} a t^2$ and the quadratic equation.
- 5) Use that time to calculate the displacement in x using $\Delta x = v_x t$

Calculations:

Setting _____, Angle _____

$v_i =$

$h =$

$t =$

$\Delta x =$

$V_y =$

$V_x =$

Setting _____, Angle _____

$v_i =$

$h =$

$t =$

$\Delta x =$

$V_y =$

$V_x =$

Setting _____, Angle _____

$v_i =$

$h =$

$t =$

$\Delta x =$

$V_y =$

$V_x =$

Setting _____, Angle _____

$v_i =$

$h =$

$t =$

$\Delta x =$

$V_y =$

$V_x =$

Results:

Setting _____, Angle _____ will hit the ground _____ from the launch point.

Setting _____, Angle _____ will hit the ground _____ from the launch point.

Setting _____, Angle _____ will hit the ground _____ from the launch point.

Setting _____, Angle _____ will hit the ground _____ from the launch point.

Day 3: Test your Predictions

Purpose: To test our predictions and catch the marbles in small cups

Procedures:

- 1) Materials: eye protection, CPO marble launcher, marble, measuring tape, small Dixie cup, tape.
- 2) Set up your marble launcher at setting 1- angle 1
- 3) Measure the predicted displacement for this arrangement and tape the Dixie cup to the floor where you think the marble will land.
- 4) Launch your marble 4 times. If after these 4 trials you haven't hit your cup, adjust the cups placement and try again.
- 5) If you had to move your cup measure its location.

Results:

Setting ____, Angle ____:

Was your prediction accurate?

If not how much did you have to adjust your cup placement?

Setting ____, Angle ____:

Was your prediction accurate?

If not how much did you have to adjust your cup placement?

Setting ____, Angle ____:

Was your prediction accurate?

If not how much did you have to adjust your cup placement?

Setting ____, Angle ____:

Was your prediction accurate?

If not how much did you have to adjust your cup placement?

Conclusion: On a separate sheet of paper write a conclusion for this lab. Be sure to discuss error.